EDUCATIONAL MOMENTS



INSPIRING CARE

How to manage patients with Corneal Staining


1. WHAT YOU NEED TO KNOW

Slit Lamp Viewing:

- 1. Parallelepiped or diffuse beam
- 2. High magnification (16 30x)

Grading Extent:


- O. None
- 1. 1–20 punctate diffuse spots
- 2. 21-40 punctate diffuse spots




- 3. Fluorescein with blue cobalt filter and yellow barrier filter
- 4. Direct illumination
- 3. >40 diffuse spots and/or coalescing patches


C. Stromal diffusion immediate but moderate

4. Dense confluent patches










Grading Depth:

A. No stromal diffusion

B. Stromal diffusion delayed (30-60 seconds)







Grade C

Grade D

Grading Position: Superior, nasal, inferior, temporal, central

Incidence:

• Common in all CL wearers (up to 60 %) but often clinically insignificant

Aetiology:

- Mechanical trauma, foreign body, damaged lens, lens edge, material stiffness, thick lens design
- Exposure disruption of tear film and subsequent desiccation
- Metabolic hypoxia, hypercapnia (tissue acidosis and desquamation of epithelial cells)
- Solution-induced corneal staining (SICS) seen 2-4h after insertion with some SiH/MPS combinations

- Some degree also seen in non CL wearers (35 %)
- Toxic care regimen hypersensitivity (1-10 % hydrogel lens wearers)
- Allergic delayed or immediate hypersensitivity
- Infectious systemic disease, poor general health (e.g. influenza, throat infection)

EDUCATIONAL MOMENTS



INSPIRING CARE

How to manage patients with Corneal Staining

1. WHAT YOU NEED TO KNOW

Symptoms:

- Can be asymptomatic depends on aetiology and severity (grade 3/4, discomfort or even pain)
- May include CL intolerance, reduced wearing time, dryness, itching
- \bullet Reduced vision if significant and over visual axis
- Lacrimation

Signs:

- Superficial punctate epithelial erosions (SPEE) (extent, depth and location depends on aetiology)
- Bulbar conjunctival hyperaemia
- Tarsal conjunctival changes
- Lacrimation

2. WHAT YOU NEED TO RECOMMEND TO YOUR PATIENTS

Recommendations:

- Manage if \geq Grade 2 or if \geq 1 grading scale increase
- Remove lenses for 24h with Grade 2, 2 3 days with Grade 3, and 7 days with Grade 4
- Consider medical intervention/treatment for Grades 3 and 4
- Ocular lubricants to reduce symptoms


- Isolate cause and manage change care system, refit/ replace lens, improve oxygen performance, blinking, rewetting drops
- For SICS: ensure a rub and rinse step (right), alter combination of SiH and MPS, switch to non-preserved solution or change to DD lens



Prognosis: Good (unless Bowman's membrane penetrated and subsequent residual scarring)

Differential diagnosis: • Infectious keratitis (right)

• Excessive eye rubbing: pressure marks



3. HOW TO FIND OUT MORE

- ► Click <u>here</u> for a refresher on slit lamp techniques
- ► Click here for the andrasko corneal staining grid
- ► Click **here** for short videos on slit lamp techniques
- ► Click **here** for further reading/references

This series is adapted from A Handbook of Contact Lens Management (3d Edition) published by THE VISION CARE INSITUTE®

MGDV : LWEV : LIPCOFV : CORNEAL MOMENT 5 : MOMENT 5 : MOMENT 7 : MOMENT 8 : MOMENT 9 : MOMENT 10 : MOMENT 11 : MOMENT 12 : MOMENT 13 : MOMENT 14 : MOMENT 15

EDUCATIONAL MOMENTS




INSPIRING CARE

How to manage patients with Corneal Staining

PATIENT CASE STUDY

When you have read this guide and our recommended resources, why not take part in THE VISION CARE INSTITUTE® self-assessment quiz to test your clinical diagnostic and management skills. Choose only one answer to each question then check the answer at the foot of the page to see whether it's correct. Good luck!



Patient GC is a 38-year-old retail assistant who has worn silicone hydrogel lenses successfully for 3 years and has been symptom-free at each aftercare visit. He comes in to the practice complaining of reduced comfort with his lenses over the past week. He recently switched from the multipurpose solution you recommended to an own-brand solution.

Questions:

- 1. What slit lamp technique would you use to examine this patient for corneal staining?
 - A. Narrow beam

C. White light

B. Low illumination

- **D.** Fluorescein and cobalt blue light
- 2. What is the best time of day to examine patients for solution-induced corneal staining?
 - A. Immediately after lens insertion

C. 6-8 hours after insertion

B. 2-4 hours after insertion

- **D.** At the end of the daily wearing time
- 3. Slit lamp examination shows this patient has approximately 30 diffuse punctate spots on the cornea. What grade would you give to this extent of staining?
 - A. Grade 1

C. Grade 3

B. Grade 2

D. Grade 4

- 4. Which of the following management options are you most likely to consider first?
 - A. Refer for medical treatment

C. Remove lenses for 24 hours then resume lens wear using your recommended solution

B. Switch to daily disposable lenses

- **D.** Ocular lubricants
- 1. The correct answer is D. Instilling fluorescein and using cobalt blue and yellow barrier filters is the best method for viewing corneal staining.
- **2.** The correct answer is B. Studies have shown that the optimum time to examine for solution-induced staining is 2-4 hours after insertion.
- **3.** The correct answer is B. A total of 21-40 punctate diffuse spots with no coalescing patches indicates Grade 2 extent of corneal staining.
- **4.** The correct answer is C. Advise the patient only to use the care product you recommend and note the brand on his contact lens specification.

